




Are you sure it is Real White??


pLUG® TA-cloning Vector Series

DCS™ Technology reduce the number of false positive colony upto 99%!
Ligation in Just 15 min!
 Short total required cloning time compatible to that of Topoisomerase-mediated cloning strategy.


Artificial T
False White



Natural
Real White!!!




pLUG® TA-cloning Vector Kit (Cat. No. 11041)



Element	Position (bp)
LacZα fragment	216 ~ 643
LacZα start codon	216
M13 reverse priming site	204 ~ 220
M13 forward priming site	439 ~ 454
Multiple cloning site	272 ~ 404
Ampicillin resistance gene	1,188 ~ 2,308
<i>fl</i> origin of replication	594 ~ 1,050
ColE1 origin of replication	2,042 ~ 2,082

pLUG®-Multi TA-cloning Vector Kit (Cat. No. 11051)



Element	Position (bp)
LacZα fragment	216 ~ 619
LacZα start codon	216
M13 reverse priming site	204 ~ 220
M13 forward priming site	371 ~ 386
Multiple cloning site	230 ~ 360
Ampicillin resistance gene	965 ~ 1,825
ColE1 origin of replication	1,826 ~ 1,866

LacZα Start

```
CAGGAAACAGCTATGACCATGATTTAGCCCAAGCTCGAARATTAACCTCACTAAAGGGAACAACAAAC
CTCCTTTGTCGATACCTGTAATAAGCGGTTTCGAGCTTTAATTGGGAGTGATTCCTCTGTTTTGCG
```

M13 reverse priming site

SacI *NotI* *XbaI* *BamHI* *SmaI* *PstI* *EcoRI*

```
TGGAGCTCCACCCGGTGGCCGCCCTCTAGAACTAGTGGATCCCGGGCTGCAGAAATTTGATA
AACCTCGAGGTGGCCCAACCGCCGGGAGATCTTGTATCACCTAGGGGGCCGACCTCTTAAGCTAT
```

HindIII *PstI* *SphI* *HindIII*

```
TCAGCTTCCAGAGCTGACCTGACCTGGCAGCTTATCGATACCTCGACCTTCAGGGGGG
AGTTCGAAAGGTCTCGTCTGCTGACCGTTCCGAATAGCTATGGCAGCTGGAAAGTCCCGCC
```

PCR Product

```
GACGTGGCAGCTTATCGATACCTCGACCTTCAGGGGGG
TCTGACCGTTCCGAATAGCTATGGCAGCTGGAAAGTCCCGCC
```

KpnI

```
CCCGGTACCCAATTCGCCCTATAGTGAAGTCTGATTAACAATTCATGGCCGCTGTTTTACACGTCGT
GGGCTATGGGTTAAGCGGGATATCACTCAGCATAATGTTAAGTACCCGGCAGCAAAATCTTCAGCA
```

M13 forward priming site

LacZα Start

```
CAGGAAACAGCTATGACCATGATTTAGCCCAAGCTCGAARATTAACCTCACTAAAGGGAACAACAAAC
CTCCTTTGTCGATACCTGTAATAAGCGGTTTCGAGCTTTAATTGGGAGTGATTCCTCTGTTTTGCG
```

M13 reverse priming site

EcoRI *SacI* *KpnI* *SmaI* *BamHI* *XbaI* *SalI*

```
CTGCAGGATCCCGGGTACCAGCTCGAATTCGAGCTTGGCACTGGCCCTCGTTTTACACGTCGT
TCTAGAGGATCCCGGGTACCAGCTCGAATTCGAGCTTGGCACTGGCCCTCGTTTTACACGTCGT
```

PCR Product

```
GACGTGGCAGCTTATCGATACCTCGACCTTCAGGGGGG
TCTGACCGTTCCGAATAGCTATGGCAGCTGGAAAGTCCCGCC
```

PstI *SphI* *HindIII*

```
TCAGCTTCCAGAGCTGACCTGACCTGGCAGCTTATCGATACCTCGACCTTCAGGGGGG
AGTTCGAAAGGTCTCGTCTGCTGACCGTTCCGAATAGCTATGGCAGCTGGAAAGTCCCGCC
```

M13 forward priming site




DCS™ Technology

The false positive colony, which is the main problem encountered in conventional TA-cloning, can be eliminated almost 99% by DCS™ Technology. **DCS™ Technology gives high percentage of true white colonies** with an anticipated recombinant plasmid.

DCS™ Technology excludes the possibility of selection of colonies having a parental plasmid used for the preparation of TA-cloning vector or a re-circulated plasmid after loss of T-overhang of TA-cloning vector.

Flow chart of TA-cloning



Component	Volume(μl)
TA-cloning vector	1
PCR product	1~4
Distilled water	Adjust to 7
5X Ligation buffer	2
T4 DNA ligase	1
Total volume	10

Characteristics

- High cloning efficiency
- High percentage of true white colony (Credible blue/white colony selection)
- Rapid procedure (Rapid ligation) : 10~15 min
- Allowing convenient sequencing
- Allowing easy re-cloning to another vector (Mirror-repeat pattern)

pLUG® TA-cloning Vector Kit (Cat. No. 11041 / 20rxn)



pLUG® TA-cloning Vector	20 μl
5X Ligation Buffer	100 μl
T4 DNA Ligase	20 μl
Distilled Water	1 ml

pLUG® - Multi TA-cloning Vector Kit (Cat. No. 11051 / 20rxn)



pLUG® - Multi TA-cloning Vector	20 μl
5X Ligation Buffer	100 μl
T4 DNA Ligase	20 μl
Distilled Water	1 ml

Distributor



iNtRON
Biotechnology

www.intronbio.com
intronbio@intronbio.com

Tel. +82-31-778-7807
Fax. +82-31-736-7245